

SAFER SIPS

STANDING TOGETHER TO STOP DRINK SPIKING



SAFER SIPS INTRODUCES CYD: CHECK YOUR DRINK.

The CYD® tests were developed and launched in 2014 as a direct response to the drink spiking problem. It has developed with expert chemists in the UK to help combat drink spiking worldwide.

WHAT DO WE TEST FOR?


CYD® tests for pharmaceutical GHB (Xyrem), Ketamine, Cocaine, Scopolamine and other amine-based drugs.

IS CYD VALIDATED?

The CYD® tests have been developed and validated by three independent bodies: Clinical Trials Laboratory Services (CTLs), University of Strathclyde and Lucideon has officially included Cocaine and Scopolamine to the Ketamine test in 2024.

Instructions:

Apply a drop of drink using a finger or a straw to the pink & yellow pads.



If either test pad changes to orange, grey or blue, **DO NOT DRINK!**



EASY TO USE

- CYD® tests are easy and quick to use with instant results.
- No Mixing. No rubbing. No waiting. Just apply a drop of drink to each pad on the test strip.
- If the colour changes to orange, grey or blue – your drink has probably been compromised or tampered with – Do Not Drink!
- Tell a trusted person, get to safety.

99% ACCURACY
97% DETECTION*

WHAT ARE THE SYMPTOMS

- Feeling unusually drunk when you haven't had any or much alcohol to feel that way.
- Sudden sweating, feeling of hot flushes or temperature change.
- Loss of feelings in your limbs.
- Loss of bowel or bladder control.
- Nausea or vomiting.
- Increased dizziness, disorientation, hallucinations
- Blurred vision.
- Memory loss or 'blackout'.

5 CHANGES TO A DRINK

- Saltiness or change in taste
- Foggy appearance
- Excessive bubbles
- Change in colour
- Sinking ice

*Some substance that are dark red or milk-based, may produce a false positive.
Bitter lemon and lime cordial may produce a false negative.